


Western Blotting: Give Your Blots a Chance to Improve


Key Things to Know Before Starting Western Blotting

Basic Knowledge

Simplified antibody structure



What is an Epitope?



Target Protein Sequence. How do I choose the right antibody for my protein?

⚠ The antibody-eliciting sequence must be present in the target protein!


Example: >sp|P83755|PSBA_ARATH Photosystem II protein D1

Signal peptide (cleaved off):


MTAILERRESESLWGRFCNWTSTENRLYGWFGVLMPTLLTATSVFIIAFAAPP
VDIGIREPVYSGSLLYGNIIISGAIPSTAAIGLHFYIWEAASVDEWLYNGGPYEL
IVLHFLGVACYMGREWLSFRLGMRPWIAVAYSAPVAAATAVFLYPIGQGSFSD
GMPLGISGTNFMVIFAENHLMHPHMLGVAGVFGGSLFSAMHGLSVTSSIRE
TTENASANEGRYRFGQEEETYNIAAHGYFGRILFYASFNNSRSLHFFLAAPVIV
GIWFTALGISTMAFNLNGFNQSVDSQGRVINTWADINRA~~NLGMVIMHERNA~~Ab to C-terminal

Antibody binding to a linear epitope requires >70% linear sequence identity

How much do you know about the protein you aim to detect?



Amount of Target Protein in a Tissue




Protein Transfer: PVDF vs. Nitrocellulose Membranes

Feature	PVDF	Nitrocellulose
Protein binding	High; hydrophobic	Moderate; hydrophobic & electrostatic
Durability	Strong; air-dried/re-wetted	Brittle & fragile; air-dried/re-wetted
Sensitivity	Ideal for low-abundance proteins	Suitable for high-abundance proteins
Pre-treatment	Methanol pre-wetting needed	No pre-wetting needed
Protein size	Broad range	Not for low MW proteins (< 20 kDa)

Secondary Antibodies


- Use at high dilutions to reduce background, cross-reactivity
- Max incubation time = 1 h/RT – do not exceed!
- Agrisera secondary Ab: high sensitivity & consistent performance



4 Key Factors for Choosing Antibodies (Recommended by Agrisera)

- Species Reactivity: Confirmed experimentally or predicted based on sequence conservation
- Validation Data: Reliable antibodies must be validated in endogenous samples for accuracy
- Tested Applications: Each technique requires a separate validation
- Publication Record: Citations on supplier websites add credibility and guide your selection

Western Blotting Workflow



Western Blot Tips & Troubleshooting

Tips


General: Get a correct band first before testing variations!

For Beginners


- Samples: WT & specific mutant
- Protein load: >50 µg/well
- Blocking: 5% nonfat milk/BSA, 1 h/RT
- Primary Ab: 1:500–1:1000, ON/4°C
- Secondary Ab: 1:25000, 1 h/RT
- Always record results!

For Advanced Users

- Review protocol and perform optimization!
- One protocol doesn't fit all antibodies



No Signal is Detected in a Total Cell Extract?




Transfer Conditions Based on Protein Size

Protein Size	Challenges	Recommended conditions			
		Membrane	MeOH	SDS	Transfer
High MW Proteins (>20 kDa)	Retained in gel; longer transfer time	0.45 µm	Low (10%)	Increased	Wet (tank)
Low MW Proteins (<20 kDa)	Easily migrates out of gel; shorter transfer time	0.20 µm	High (20%)	Decreased	Semi-dry

MeOH is toxic! EtOH is a safer, equally effective alternative to MeOH

Perform wet and semi-dry transfers under cold conditions to prevent heat buildup

Can Primary Antibodies be Reused?



How to Determine Gel Loading?


Loading Method	Recommendation	Notes
µg of Protein	Recommended	Best for accurate loading
µg of Chlorophyll	Recommended	Suitable for plant samples
Volume (µl)	Not recommended	Can lead to inconsistent results

For example, there are reagents which allow fast protein quantification and are compatible with detergents, salts, and chelators present in the extraction buffer

For photosynthetic tissues


Western Blotting: Give Your Blots a Chance to Improve. For more detailed information, scan QR codes and check specific references and methods therein. Send questions, comments, and suggestions to Joanna Porankiewicz-Asplund (joanna@agrisera.com). **Abbreviations:** AA, amino acids; Ab, antibody; BCA, bicinchoninic acid; BSA, bovine serum albumin; ChIP, chromatin immunoprecipitation; EtOH, ethanol; HRP, horseradish peroxidase; IgG, immunoglobulin G; IHC, immunohistochemistry; IP, immunoprecipitation; KO/KD, knockout/knockdown; MeOH, methanol; MW, molecular weight; ON, overnight; PAXDb, protein abundance database; PBST, phosphate-buffered saline with Tween® 20 detergent; PSI, photosystem I; ppm, parts per million; PVDF, polyvinylidene fluoride; RT, room temperature; SDS-PAGE, sodium dodecyl sulfate-polyacrylamide gel electrophoresis; TBST, tris-buffered saline with Tween® 20 detergent; WB, Western Blot; WT, wild type. **Notes:** Protein complexes were generated using Protein Imager using coordinates deposited at PDB with IDs 3jcu and 7ou1. **Acknowledgements:** We thank Anna Węgrzyn, David Alabadi, Dhruv Agrawal, Edith Kalén, Louice Lindroth, Mariana Barber, Mika Teranishi, Pushan Bag, Sebastian Hoernstein, and Ximena Anleu Gil for valuable comments and corrections. We are highly grateful to Agrisera for sponsoring poster design, printing, and free distribution at conferences around the world. Citation: Porankiewicz-Asplund J, Shevela D, et al (2025) Western Blotting: Give Your Blots a Chance to Improve, Agrisera Educational Poster 9. doi:10.6084/m9.figshare.2908685.

Quick Protein Transfer Quality Control



Is Sample Denaturation Above 90°C a Golden Standard?

Answer: Not always! High temperatures can cause protein aggregation, and signal loss in Western Blot, while membrane proteins may never enter the gel



Antibody Validation Checklist


Genetic Validation (Direct approach)

- KO/KD Control: Target protein signal disappears or decreases in KO/KD samples
- IP-MS: Confirms antibody binds the correct protein

Biochemical & Experimental Validation (Indirect approach)

- Correct MW: Band appears at expected MW, considering modifications
- Low Cross-Reactivity: Minimal non-specific binding
- Cellular Localization: Signal detected in the correct cellular fraction
- Expression-Based Validation: Signal changes with known regulatory conditions
- Independent Antibody Confirmation: Another validated antibody detects the same band
- Pre-immune Serum Control: No/minimal signal in pre-immune serum
- Peptide Neutralization: Signal reduces with target peptide pre-incubation

How to Cite an Antibody? Why it Matters?



How to Get a Strong Signal from a Weak Antibody?

