This product is **for research use only** (not for diagnostic or therapeutic use)

contact: support@agrisera.com

Agrisera AB | Box 57 | SE-91121 Vännäs | Sweden | +46 (0)935 33 000 | www.agrisera.com

Product no AS20 4394

Aflatoxin M1 | ELISA quantitation kit

Background

This ELISA assay utilizes the principle of competitive binding to measure the concentration of Aflatoxin M1 in analyzed samples. The anti-Aflatoxin M1 antibody has been pre-coated on the surface of the reaction wells. Samples, containing an unknown amount of Aflatoxin M1, or standards are added to the appropriate plate wells together with Horseradish Peroxidase (HRP) conjugated Aflatoxin M1. During incubation, the competitive inhibition reaction occurs between (HRP) conjugated Aflatoxin M1 and Aflatoxin M1 in standards and samples. Substrate solution is added to the wells and the color develops in opposite to the amount of Aflatoxin M1 in the sample or standards. Reaction is stopped and the intensity of the color is measured at 450 nm.

The Aflatoxin M1 ELISA kit reagents

Assay plate	1 (96 wells)
Standard	5 x 1.5 ml
HRP-conjugate	1 x 15 ml
Substrate A	1 x 7 ml
Substrate B	1 x 7 ml
Stop Solution	1 x 7 ml
Wash Buffer (10x concentrate)	1 x 40 ml
Sample Diluent 1	1 x 10 ml
Sample Diluent 2	1 x 20 ml

Sensitivity: ≥ 5 ppt **Detection range:** 5-135 ppt

Recovery rate: cheese, cream, milk, milk powder, raw milk, yoghurt: 95% ± 30%

Limit of detection: cheese: 15 ppt, cream: 45 ppt, raw milk, milk: 5 ppt, milk powder: 25 ppt, yoghurt: 10 ppt

Detection wavelength: 450 nm

Intra-assay precision (within assay): CV%<8% Intra-assay precision (between assays): CV%<10%

Storage: 2-8°C

Cross-reactivity: Aflatoxin M1 (AFM1) 100%, Aflatoxin B1 (AFB1) < 1% **Sample type:** cheese, cream, milk, milk powder, raw milk, yoghurt

This product is **for research use only** (not for diagnostic or therapeutic use)

Agrisera AB | Box 57 | SE-91121 Vännäs |Sweden | +46 (0)935 33 000 | www.agrisera.com

Not provided in the kit but required to conduct the test

Microplate reader capable of measuring absorbance at 450 nm, with the correction wavelength set at 540 nm or 570 nm | Incubator with stable incubation temperature up to 25°C | Squirt bottle, manifold dispenser or automated microplate washer | Rotary evaporator or nitrogen gas | Analytical balance, 2 decimal place | Single-channel micropipette (20-200 μ l, 100-1000 μ l) | Multi-channel micropipette (30-300 μ l) | Methanol | Ethanol

contact: support@agrisera.com

Important notes

The kit should not be used beyond the expiration date on the kit label. Reagents from different lots should not be mixed. In case analyzed samples are generating values which are higher than the highest standard, samples should be diluted, and the assay repeated. The assay may vary due to factors such as operator, pipetting technique, washing technique, incubation time, temperature, and kit expiry date. This assay is designed to eliminate the interference caused by soluble receptors, binding proteins, and other components present in biological samples. Until all factors have been tested in the Immunoassay, the possibility of interference cannot be excluded.

Samples: Samples should be added very carefully to avoid foaming and touching the well wall. For each step in the procedure, the total dispensing time for the addition of reagents or samples to the assay plate should not exceed 10 min. Cross contamination should be avoided by changing pipette tips between additions of each standard, sample and reagent.

Incubation: Proper adhesion of plate sealers during each incubation step has to be applied to ensure accurate results. Strips should not be let dry at any time. Incubation temperature and time must be kept constant with the optimal temperature being 25°C.

Washing: This step is critical, and the complete removal of liquid at each step is essential for a good performance of this assay. Any drops and fingerprints must be removed from a bottom of the plate. Insufficient washing will result in poor precision. In case of automated washing, addition of a 30 s soak period, following the addition of wash buffer and plate rotation by 180 degrees between each wash step, may improve assay precision.

Reaction development: Substrate should change from colorless or light blue to gradations of blue. In case the color is too intense Stop Solution should be added. Substrate solution can be easily contaminated. The liquid should remain colorless or light blue until added to the plate. Keep this solution away from any light.

Stop solution: The reaction is developed correctly if color development in the wells is from blue to yellow upon addition of the Stop Solution. Green color indicates that the Stop Solution was not thoroughly mixed with the substrate. In case of absorbance value of standard solution 0 being less than 0.5 this indicates its degeneration.

Safety measures

The Aflatoxin M1 ELISA quantitation kit Stop Solution is an acidic solution. Please use suitable protective gear for your eyes, hands, face, and clothes when handling this item.

This product is **for research use only** (not for diagnostic or therapeutic use)

Agrisera AB | Box 57 | SE-91121 Vännäs | Sweden | +46 (0)935 33 000 | www.agrisera.com

Procedure instructions

Please read the whole manual carefully before proceeding with your experiment.

Before starting

As a recommendation, all samples and standards should be made in duplicates for this assay. Do not exceed handling time of 10 minutes per step. Always use a plate cover during all steps. Do not let the strips dry at any time during the assay.

contact: support@agrisera.com

Sample preparations

Raw milk/milk

Use 100 µl of homogenized sample for the analysis. Dilution factor of the sample: 1.

Milk powder

Add 4 ml of water to 1.0 ± 0.1 g of homogenized sample and vortex until the sample is completely dissolved. Heat in a water bath at 60°C for 5 min. Centrifuge at a minimum of 4000 rpm for 5 min. Use 100 μ l of the under layer for further analysis. Dilution factor of the sample: 5.

Yoghurt

Add 0.1 ml of Sample Diluent 1 to 1.0 \pm 0.1 g of homogenized sample. Add 1 ml of deionized/distilled water and vortex for 3 min. Centrifuge at a minimum of 4000 rpm for 5 min. Use 100 μ l of the under layer for further analysis. Dilution factor of the sample: 2.

Cream

Add 2 ml of Ethanol (30%) to 1.0 ± 0.1 g of homogenized sample and vortex for 3 min. Heat in a water bath at 60°C for 5-10 min, until completely dissolved. Centrifuge at a minimum of 4000 rpm for 5 min. Take 100 μ l of the middle layer and add 200 μ l Sample Diluent 2. Vortex well. Use 100 μ l for further analysis. Dilution factor for the sample: 9.

Cheese

Add 2 ml of Methanol (50%) to 1.0 ± 0.1 g of homogenized sample. Heat in a water bath at 60°C for 5-10 min, until completely dissolved. Vortex for 3 min. Centrifuge at a minimum of 4000 rpm for 5 min. Take 500 μ l of the middle layer at 60°C by nitrogen or rotary evaporator. Add 500 μ l of the Sample Diluent 2 and vortex well. Use 100 μ l for further analysis. Dilution factor for the sample: 3.

This product is **for research use only** (not for diagnostic or therapeutic use)

contact: support@agrisera.com

Agrisera AB | Box 57 | SE-91121 Vännäs | Sweden | +46 (0)935 33 000 | www.agrisera.com

Manual

- 1. Bring all reagents and samples to room temperature before use (keep in RT for at least 30 min).
- 2. Centrifuge the samples after thawing.
- 3. Prepare the reagents and samples.
 - a) Prepare the Washing Buffer (1x). In case crystals have formed, warm the concentrate to room temperature and mix gently until the crystals have dissolved. Prepare 90 ml of Wash Buffer (1x) by diluting 10 ml of Washing Buffer Concentrate (10x) with deionized or distilled water.
 - b) Prepare the Standards. Centrifuge the standard vials at 6000-10000 rpm for 30 s. The most concentrated standard, S4, is at 135 ppt, while the zero standard, S0, is at 0 ppt.

Tube	S4	S3	S2	S1	S0
ppt	135	45	15	5	0

- 4. Determine the number of wells to be used. Unused wells should be put back into the Ziplock pouch and stored at 4°C.
- 5. Set one blank well without any solution. Add 100 μ l of standard or sample to each well. Add samples gently to avoid foaming and be careful not to touch the well walls and incubate at 25°C for 30 min.
- 6. Aspirate the wells and wash the plate by filling the wells with Wash Buffer (1x, 250 μl in each well) using a squirt bottle, multi-channel pipette or auto-washer. Let it stand for 15-30 s, remove the liquid, and repeat the washing process four times for a total of five washes. After the fifth wash, remove all liquid and invert the plate against clean paper towels. Remove water and fingerprint on the bottom of the plate to avoid falsely reading results. The wash procedure is critical. Complete removal of liquid at each step is essential.
- 7. Add 100 μ l of HRP-conjugate (1x) to each well (except the blank well). Mix well and incubate at 25°C for 15 min.
- 8. Repeat the washing procedure for five times as in step 6.
- 9. Add 50 μ l of Substrate A and 50 μ l of Substrate B to each well and mix well. Incubate at 25°C for 15 min in the dark. (Check the color development from colorless or light blue to blue every 10 min, stop in advance if the color is too deep.)
- 10. Add 50 μ l of Stop Solution to each well in the same order as the substrates and mix thoroughly. The color will change from blue to yellow.
- 11. Read the optical density of each well within 5 minutes, at 450 nm.

Sample	ppt	Average
S4	135	1.653
S3	45	1.385
S2	15	1.056
S1	5	0.826
S0	0	0.453

This product is **for research use only** (not for diagnostic or therapeutic use)

contact: support@agrisera.com

Agrisera AB | Box 57 | SE-91121 Vännäs | Sweden | +46 (0)935 33 000 | www.agrisera.com

Results

There are two methods to judge the results: the first one (A) is the rough judgment, while the second (B) is the quantitative determination.

Note that the OD value of the sample has a negative correlation with Aflatoxin M1 in the sample.

A: Compare the sample average absorbance values with standards values, the Aflatoxin M1 concentration in the samples can be concluded. For example, the absorbance value of sample 1 is 0.5, the absorbance value of sample 2 is 1.0; absorbance values of standard are: 1.653, 1.385, 1.056, 0.826, 0.453 and the corresponding concentrations are: 0 ppt, 5 ppt, 15 ppt, 45 ppt, 135 ppt, then the Aflatoxin M1 in sample 1 and sample 2 are 45 ppt-135 ppt and 15 ppt-45 ppt. Lastly, the reader is multiplied by the corresponding dilution factor of each sample and the actual concentration of sample is obtained.

B: The software offered together will facilitate the calculation process, it is suitable for accurate and fast analysis of large-scale samples. Please contact us for more information. Note: Discard substrate with any color that indicates the degeneration of this solution. When the absorbance value of standard solution 0 is less than 0.5, this indicates its degeneration. The optimum reaction temperature is 25°C. A temperature that is too high or too low will result in changes in the absorbance value and detecting sensitivity.

Related products

Aflatoxin M1 (50 ug)

AS14 2801 | Reactivity: Aflatoxin M1

Aflatoxin B1 | ELISA quantitation kit

AS20 4397 | Reactivity: Aflatoxin B1 from edible oil, feedstuff, feed, grain