

This product is for research use only (not for diagnostic or therapeutic use)

contact: support@agrisera.com

Agrisera AB | Box 57 | SE-91121 Vännäs | Sweden | +46 (0)935 33 000 | www.agrisera.com

Product no AS15 2843

Llama IgG - AffinityPrep Resin

Product information

Host Llama

Purity Purified Llama IgG coupled to 4% agarose beads in 10 mM Sodium Phosphate, 0.15 M Sodium Chloride, pH 7.2.

Contains 0.05 % sodium azide.

Format Liquid

Quantity 2 ml

Storage 2-8 °C, this product should not be frozen. The shelf lite is 1 year from date of receipt. Prepare working dilution prior to

use and then discard.

Additional information 10 mM Sodium Phosphate, 0.15 M Sodium Chloride, pH 7.2, 0.05% (w/v) Sodium Azide.

This resin is suitable for removal of proteins recognizing Llama IgG from liquid solutions, including traditional

Immunoabsorptions and Immunoprecipitation protocols.

8-10 mg protein/ml resin

Application information