

This product is for research use only (not for diagnostic or therapeutic use)

contact: support@agrisera.com

Agrisera AB | Box 57 | SE-91121 Vännäs | Sweden | +46 (0)935 33 000 | www.agrisera.com

Product no AS18 AgriCtBlock-100

AgriCoat-Block | Non mammalian Protein Coat Block (5X) (100 ml)

Product information

Agrisera Non Mammalian Protein Coat Block is designed to stabilize coated ELISA plates and other solid phase substrates. The reagent was developed to substantially reduce non specific antibody binding and eliminate non specific protein interactions. It provides structural suport and crowding agents and minimizes background interactions and issues.

Format Liquid

Quantity 100 ml

Recommendation for use:

Additional information

Can be used for antigens and antibody coating, blocking or as a diluent. Blocking and stabilization in one step may allow to eliminated an additional blocking step.

Need to be equlibrated to ambient temperature before use.

To prepare a working solution: dilute with sterile, deionized water.

This product synergizes well with AS15 HRPstab-100 Agrisera HRP-Conjugate/Protein Stabilizer (100 ml)

Application information

Recommended dilution

Provided as a 5x concentrate when diluted 10x, this reagent can serve as an excellent general, neutral buffered dilluent

Stock solution is viscous. Equilibrate to room temperature before use.

Additional information

For stabilization of coated ELISA plates and other solid phase substrates.

Contains anti-micorbial agent.