

This product is for research use only (not for diagnostic or therapeutic use)

contact: support@agrisera.com

Agrisera AB | Box 57 | SE-91121 Vännäs | Sweden | +46 (0)935 33 000 | www.agrisera.com

Product no AS15 TMB-HRP-10 AgriseraTMB HRP Substrate (10 ml)

Product information

Quantity	10 ml
Storage	Stable for 3 years at +2-8° C or at room temperature up to six months, and can also be frozen. Storge at room temperature for 1 month results in ca, 2 % activity loss. Protect the substrate from exposure to sunlight.
Additional information	One component ready to use reagent. It contains 2.5 mMol L-1 TMB and Hydrogen Peroxide in a buffer at pH 3.6. It also contains non-toxic stabilizers.
	How to work with this reagent?
	 Pour out amount of reagent which is needed for your assay. Do not our back any unused substrate back to the original bottle. Any trace metals, HRP or other contaminant may be reason for the substrate to turn blue.

Application information