

This product is for research use only (not for diagnostic or therapeutic use)

contact: support@agrisera.com

Agrisera AB | Box 57 | SE-91121 Vännäs | Sweden | +46 (0)935 33 000 | www.agrisera.com

Product no AS18 AgriCtBlock-1L AgriCoat-Block | Non mammalian Protein Coat Block (5X) (1L ml)

Product information

Agrisera Non Mammalian Protein Coat Block is designed to stabilize coated ELISA plates and other solid phase substrates. The reagent was developed to substantially reduce non specific antibody binding and eliminate non specific protein interactions. It provides structural suport and crowding agents and minimizes background interactions and issues.

Format	Liquid
Quantity	1L
Additional information	Recommendation for use:
	Can be used for antigens and antibody coating, blocking or as a diluent. Blocking and stabilization in one step may allow to eliminated an additional blocking step. Need to be equilbrated to ambient temperature before use. To prepare a working solution: dilute with sterile, deionized water. This product synergizes well with <u>AS15 HRPstab-100</u> Agrisera HRP-Conjugate/Protein Stabilizer (100 ml)
Application information	
Performended dilution	Provided as a Ex concentrate when diluted 10x, this regress can carve as an excellent general neutral buffered diluce

 Recommended dilution
 Provided as a 5x concentrate when diluted 10x, this reagent can serve as an excellent general, neutral buffered dilluent

 Additional information
 Stock solution is viscous. Equilibrate to room temperature before use.

 For stabilization of coated ELISA plates and other solid phase substrates.
 Contains anti-micorbial agent.