

This product is for research use only (not for diagnostic or therapeutic use)

contact: support@agrisera.com

Agrisera AB | Box 57 | SE-91121 Vännäs | Sweden | +46 (0)935 33 000 | www.agrisera.com

# Product no AS18 SecondaryECL

## Agrisera Secondary/ECL offer

### **Product information**

Offer contains the following products:

#### 1. Secondary antibody

Goat anti-Rabbit IgG, HRP-conjugated-trial sample(10 µl), to be used at dilutions from 1: 25 000 to 1: 50 000, 1h/RT incubation. Tested on a wide range of species including: dicots, monocots, moss, algae, cyanobacteria and diatoms. High dilution of this antibody gives the benefit of low background signal.

Enough for 50 Western blots (10 ml incubation volume/test).



### Agrisera ECL kit (Bright/SuperBright) trial pack

ECLBright - low pico to mid femtogram detectio, for detection of protein of moderate and high expression ECLSuperBright - extreme low femtogram detection, for detection of protein with low expression level



This set combines Agrisera primary antibody with matching effective secondary antibody and high and medium sensitivity chemiluminescent detection reagents.

# **Application information**