

This product is for research use only (not for diagnostic or therapeutic use)

contact: support@agrisera.com

Agrisera AB | Box 57 | SE-91121 Vännäs | Sweden | +46 (0)935 33 000 | www.agrisera.com

Product no AS15 TMB-HRP-10

## AgriseraTMB HRP Substrate (10 ml)

## **Product information**

Quantity 10 ml Storage

Stable for 3 years at +2-8° C or at room temperature up to six months, and can also be frozen. Storge at room temperature for 1 month results in ca, 2 % activity loss. Protect the substrate from exposure to sunlight.

**Additional information** 

One component ready to use reagent. It contains 2.5 mMol L-1 TMB and Hydrogen Peroxide in a buffer at pH 3.6. It also contains non-toxic stabilizers.

## How to work with this reagent?

- Pour out amount of reagent which is needed for your assay.
- Do not our back any unused substrate back to the original bottle.
- Any trace metals, HRP or other contaminant may be reason for the substrate to turn blue.

## **Application information**