

This product is for research use only (not for diagnostic or therapeutic use)

contact: support@agrisera.com

Agrisera AB | Box 57 | SE-91121 Vännäs | Sweden | +46 (0)935 33 000 | www.agrisera.com

Product no AS08 281

95 kDa Lcm (phycobilisome - to- thylakoid core linker with phycocyanobilin chromophore) Product information

Immunogen native protein purified from phycobilisomes of *Porphyridium cruentum*

Host Rabbit

Clonality Polyclonal

Purity Serum

Format Lyophilized

Quantity 200 μl

Reconstitution For reconstitution add 200 μl of sterile water

Storage Store lyophilized/reconstituted at -20 °C; once reconstituted make aliquots to avoid repeated freeze-thaw cycles. Please remember to spin the tubes briefly prior to opening them to avoid any losses that might occur from material adhering to

the cap or sides of the tube.

Additional information Anabaena sp. sample can serve as a negative control, as under nitrogen deficient conditions phycobiliproteins are going to be lost.

Overall sample quality is of crucial importance and in older or not properly stored samples, phycobiliproteins will

undergo proteolytic degradation.

Application information

Recommended dilution 1:300-1:1000 (WB)

Expected | apparent 95 kDa (*Porphyridium cruentum*), 75-120 kDa (other species)

Confirmed reactivity Anabaena variablis, Aphanathece halophytica, Fremyella diplosiphon (Microchaete diplosiphon CCALA 811),

Mastigocladus laminosus, Nostoc sp., Synechococcus sp. PCC 7002, ATCC 27264, PR6), Synechococcus elongatus (strain PCC 7942), Anacystis nidulans R2), Synechococcus leopolensis (..elangatus PCC 6301), Synechecoccus lividus, Spirulina platensis (Athrospira platensis), Tolypothrix tenuis; red alga: Porphyridium cruentum (purpureum), P.

aeruginuso, Rhodosorus marinus, Griffithsia monilis

Predicted reactivity Algae (red), Cyanobacteria

Species of your interest not listed? Contact us

Not reactive in No confirmed exceptions from predicted reactivity are currently known

Selected references Redlinger & Gantt (1982) A Mr 95,000 polypeptide in Porphyridium cruentum phycobilisomes and thylakoids: PNAS

79:5542

Redlinger & Gantt (1981). Phycobilisome structure of Porphyridium cruentum. Plant Physiol. 68:1375.