

This product is for research use only (not for diagnostic or therapeutic use)

contact: support@agrisera.com

Agrisera AB | Box 57 | SE-91121 Vännäs | Sweden | +46 (0)935 33 000 | www.agrisera.com

Product no AS06 179

CDJ1 | chloroplast DnaJ homolog 1

Product information

Immunogen fusion protein encoding region for the *Chlamydomonas rerinhardtii* N-terminal part of CDJ1 Q66YD3, overexpressed in

Host Rabbit

Clonality Polyclonal

Purity Serum

Format Lyophilized

Quantity 100 μl

Reconstitution For reconstitution add 100 μl of sterile water

Storage Store lyophilized/reconstituted at -20°C; once reconstituted make aliquots to avoid repeated freeze-thaw cycles. Please remember to spin the tubes briefly prior to opening them to avoid any losses that might occur from material adhering to

the cap or sides of the tube.

Application information

Recommended dilution 1, 3000 (WB)

Expected | apparent MW 40 | 40 kDa

Confirmed reactivity | Chlamydomonas reinhardtii

Predicted reactivity Ostreococcus sp.

Species of your interest not listed? Contact us

Not reactive in Higher plants

Additional information Western blot has been done on 10-50 µg load per well; immunoprecipitation has been done on 5-50 µl

Selected references Willmund et al (2008). The chloroplast DnaJ homolog CDJ1 of Chlamydomonas reinhardtii is part of a multichaperone

complex containing HSP70B, CGE1, and HSP90C. Plant Physiol. 148: 2070-2082.