## 

This product is for research use only (not for diagnostic or therapeutic use)

contact: support@agrisera.com

Agrisera AB | Box 57 | SE-91121 Vännäs | Sweden | +46 (0)935 33 000 | www.agrisera.com

## Product no AS10 938 COX11 | Cytochrome c oxidase assembly protein

## **Product information**

ImmunogenRecombinant fragment of Chlamydomonas reinhardtii COX11 UniProt: A8JCS0HostRabbitClonalityPolyclonalPuritySerumFormatLyophilizedQuantity200 μlReconstitutionFor reconstitution add 200 μl of sterile waterStorageStore lyophilized/reconstituted at -20°C; once reconstituted make aliquots to avoid repeated freeze-thaw cycles. Please<br/>remember to spin the tubes briefly prior to opening them to avoid any losses that might occur from material adhering to<br/>the cap or sides of the tube.

## **Application information**

| Recommended dilution      | 1 : 1000 (WB)  |
|---------------------------|--|
| Expected   apparent<br>MW | 33,7 kDa   |
| -                         | Chlamydomonas reinhardii   |
| Predicted reactivity      | Physcomitrium patens<br>Species of your interest not listed? <u>Contact us</u> |
| Not reactive in           | Arabidopsis thaliana   |
| Additional information    | Antibody works on purified, Chlamydomonas reinhardtii mitochondria.            |