Agrisera Part of Olink[®] Group \$

This product is for research use only (not for diagnostic or therapeutic use)

Agrisera AB | Box 57 | SE-91121 Vännäs | Sweden | +46 (0)935 33 000 | www.agrisera.com

Product no AS15 3073

Anti-DnaK

Product information

Immunogen	purified DnaK of Yersinia pseudotuberculosis YPIII/pIB102
Host	Rabbit
Clonality	Polyclonal
Purity	Serum
Format	Lyophilized
Quantity	50 μl
Reconstitution	For reconstitution add 50 µl of sterile water
Storage	Store lyophilized/reconstituted at -20 °C; once reconstituted make aliquots to avoid repeated freeze-thaw cycles. Please remember to spin the tubes briefly prior to opening them to avoid any losses that might occur from material adhering to the cap or sides of the tube.

Application information

1 : 1000 (WB)
Yersinia pseudotuberculosis YPIII/pIB102
Chlamydomonas reinhardtii
Species of your interest not listed? Contact us

Not reactive in No confirmed exceptions from predicted reactivity are currently known

contact: support@agrisera.com