

This product is for research use only (not for diagnostic or therapeutic use)

contact: support@agrisera.com

Agrisera AB | Box 57 | SE-91121 Vännäs | Sweden | +46 (0)935 33 000 | www.agrisera.com

Product no AS99 001-ALP

Anti-KLH | Keyhole limpet hemocyanin, ALP-conjugated (40 µg)

Product information

Immunogen Purified keyhole limpet hemocyanin (<u>KLH</u>), whole molecule.

Host Rabbit

Clonality Polyclonal

Purity Immunogen affinity purified serum in PBS pH 7.4, conjugated to ALP.

Format Liquid

Quantity 40 μg

Storage Storage Store at 4°C for 12-18 months. A preservative may be added for long time storage up to 2 years.

Additional information Protein present in plant vascular tissue (xylem and vascular cambium) is detected by anti-KLH antibodies (Höglund et

al. 2002) which might lead to false results in IL when using anti-peptide antibodies generated to KLH-conjugated peptide. Further information about it can be found here.

peptide. I diffier information about it carri

Application information

Recommended dilution 1:5 000 (ELISA), 1:1000 (IL), 1:5 000 (WB)

Expected | apparent

ca. 400 kDa/subunit

Confirmed reactivity Megathura crenulata - most commonly used carrier protein

Not reactive in No confirmed exceptions from predicted reactivity are currently known

Additional information Antibody can be used as a negative control to determine if observed signal is generated by anti-KLH or anti-peptide

antibodies. Due to its large size KLH protein will be very difficult to separate on SDS-PAGE.

Optimal working dilution has to be determined by end user.