This product is for research use only (not for diagnostic or therapeutic use)

Agrisera AB | Box 57 | SE-91121 Vännäs | Sweden | +46 (0)935 33 000 | www.agrisera.com

Product no AS10 1496 Rabbit anti-Guinea pig IgG (H&L), HRP conjugated

Product information

| Immunogen | Purified guinea pig IgG (H&L) <u>NP_001166520</u> |
|-------------------------|---|
| Host | Rabbit |
| Clonality | Polyclonal |
| Purity | Immunogen affinity purified rabbit IgG. |
| Format | Lyophilized |
| Quantity | 1 mg |
| Reconstitution | For reconstitution add 1,1 ml of sterile water, Let it stand 30 minutes at room temperature to dissolve, Prepare fresh working dilutions daily |
| Storage | Store lyophilized material at 2-8°C. For long time storage after reconstitution, dilute the antibody solution with glycerol to a final concentration of 50% glycerol and store as liquid at -20°C, to prevent loss of enzymatic activity. For example, if you have reconstituted 1 mg of antibody in 1,1 ml of sterile water add 1,1 ml of glycerol. Such solution will not freeze in -20°C. If you are using a 1:5000 dilution prior to diluting with glycerol, then you would need to use a 1:2500 dilution after adding glycerol. Prepare working dilution prior to use and then discard, Be sure to mix well but without foaming. |
| Additional information | Purity of this preparation is > 95% based on SDS-PAGE. Antibody concentration is 1.0 mg/ml. Antibody is supplied in 10 mM sodium phosphate, 0.15 M sodium chloride, pH 7.2.1 % (w/v) B, Protease/IgG free. Contains 0.1 % (v/v) Kathon CG as preservative of bacterial growth. |
| | Based on immunoelectrophoresis, this antibody reacts with: heavy chains on guinea pig IgG, light chains on all guinea pig immunoglobulins. Based on immunoelectrophoresis, no reactivity is observed to: non-immunoglobulin guinea pig serum proteins |
| Application information | |

Application information

Recommended dilution This conjugate is suitable for all immunoassay applications. The optimal working dilution should be determined by the investigator,

contact: support@agrisera.com

