

This product is for research use only (not for diagnostic or therapeutic use)

contact: support@agrisera.com

Agrisera AB | Box 57 | SE-91121 Vännäs | Sweden | +46 (0)935 33 000 | www.agrisera.com

Product no AS15 2997

Anti-SPSC | Sucrose phosphate synthase isoform C

Product information

KLH-conjugated synthetic peptide derived from Arabidopsis thaliana SPS4 sequence, Uniprot: F4JLK2, TAIR: Immunogen

Host Rabbit

Clonality Polyclonal

Purity Immunogen affinity purified serum in PBS pH 7.4.

Format Lyophilized

Quantity 50 μg

Reconstitution For reconstitution add 50 μl of sterile water

Storage | Store lyophilized/reconstituted at -20°C; once reconstituted make aliquots to avoid repeated freeze-thaw cycles. Please remember to spin the tubes briefly prior to opening them to avoid any losses that might occur from material adhering to

the cap or sides of the tube.

Application information

Recommended dilution 1: 1000 (WB)

Expected | apparent

118.88 kDa MW

Predicted reactivity Corchorus olitorius (A0A1R3I4K1), Gossypium hirsutum (A0A1U8K702)

Species of your interest not listed? Contact us

Not reactive in No confirmed exceptions from predicted reactivity are currently known

Application example



Proteins were isolated from Arabidopsis thaliana wilde-type and SPS mutants were solubilized with 3X LB (6 M urea, 12% SDS, 30% glycerol, 100 mM DTT, 150 mM Tris pH7.0, 0.8% Comassie G-250). 10 µg of proteins were loaded into each lane and separated on 10% SDS-PAGE, and then blotted overnight onto PVDF membrane. Membranes were blocked with milk powder for 2 h and then incubated in the primary antibody solution at 1: 1 000 overnight at 4°C with agitation, which was then decanted and the membrane was washed 3 times for 5 min in TBST. Membrane was incubated at RT for 1 hour in 1:10 000 Agrisera goat anti-Rabbit secondary antibody (AS09 602), followed by washing steps as above. Membrane was developed for 2 min with chemiluminescent detection reagent, according to the manufacturer's instructions and recorded using FujiFilm CCD camera with 10 s increment time for around 5 min.