

Product no **AS22 4730****Anti-PsaD | PSI-D subunit of photosystem I (algal)****Product information**

Anti-PsaD | PSI-D subunit of photosystem I (algal)

Immunogen	PsaD of <i>Chlamydomonas reinhardtii</i> , UniProt: Q39615
Host	Rabbit
Clonality	Polyclonal
Purity	Serum
Format	Lyophilized
Quantity	50 µl
Reconstitution	For reconstitution add 50 µl of sterile water
Storage	Store lyophilized/reconstituted at -20°C; once reconstituted make aliquots to avoid repeated freeze-thaw cycles, Please, remember to spin tubes briefly prior to opening them to avoid any losses that might occur from lyophilized material adhering to the cap or sides of the tubes

Application information

Recommended dilution	1: 5000 - 1 : 10 000 (WB)
Expected apparent MW	21.3 kDa
Confirmed reactivity	<i>Chlamydomonas reinhardtii</i>
Additional information	Can be provided with ProClin, if requested