

Product no **AS05 056****ExoS | Exoenzyme S****Product information**

<b>Immunogen</b>	amino acids 366 to 453 of <a href="#">PA3841</a> of ADP-ribosylating enzyme - Exoenzyme S overexpressed in a GST fusion. Afterwards cleaved with a help of trombin and separated on a polyacrylamide gel. Gel piece has been used for immunizations.
<b>Host</b>	Chicken
<b>Clonality</b>	Polyclonal
<b>Purity</b>	Purified, total IgY (chicken egg yolk immunoglobulin) in PBS pH 8. Contains 0.02 % sodium azide.
<b>Format</b>	Liquid
<b>Quantity</b>	100 µl
<b>Storage</b>	Store at 4°C; make aliquots to avoid working with a stock. Please remember to spin the tubes briefly prior to opening them to avoid any losses that might occur from material adhering to the cap or sides of the tube.

**Application information**

<b>Recommended dilution</b>	1 : 5000 (WB)
<b>Expected   apparent MW</b>	48 kDa
<b>Confirmed reactivity</b>	<i>Pseudomonas aeruginosa</i>
<b>Predicted reactivity</b>	<i>Pseudomonas aeruginosa</i>
<b>Not reactive in</b>	No confirmed exceptions from predicted reactivity are currently known
<b>Selected references</b>	<p><a href="#">Feng et al. (2019): Tanshinones: First-in-Class Inhibitors of the Biogenesis of the Type 3 Secretion System Needle of Pseudomonas aeruginosa for Antibiotic Therapy. ACS Cent. Sci.2019.</a></p> <p><a href="#">Anantharajah et al. (2017). Salicylidene acylhydrazides and hydroxyquinolines act as inhibitors of type three secretion systems in Pseudomonas aeruginosa by distinct mechanisms. Antimicrob Agents Chemother. 2017 Apr 10. pii: AAC.02566-16. doi: 10.1128/AAC.02566-16.</a></p> <p><a href="#">Anantharajah et al. (2016). Inhibition of the Injectisome and Flagellar Type III Secretion Systems by INP1855 Impairs Pseudomonas aeruginosa Pathogenicity and Inflammasome Activation. J Infect Dis. 2016 Jul 13. pii: jiw295.</a></p>