Aa sera Part of Olink® Group

This product is for research use only (not for diagnostic or therapeutic use)

contact: support@agrisera.com

Agrisera AB | Box 57 | SE-91121 Vännäs | Sweden | +46 (0)935 33 000 | www.agrisera.com

Product no AS14 2799 Anti-AGO3-PAZ | Argonaute 3 PAZ domain (Chlamydomonas)

Product information

Immunogen KLH-conjugated peptide derived from C-terminal of AGO3 protein of Chlamydomonas reinhardtii Host Rabbit Clonality Polyclonal **Purity** Immunogen affinity purified serum in PBS pH 7.4. Format Lyophilized Quantity 50 µg Reconstitution For reconstitution add 50 µl of sterile water Storage Store lyophilized/reconstituted at -20°C; once reconstituted make aliquots to avoid repeated freeze-thaw cycles. Please remember to spin the tubes briefly prior to opening them to avoid any losses that might occur from material adhering to the cap or sides of the tube.


Application information

Recommended dilution 1: 2500 (WB) Expected | apparent MW Confirmed reactivity Chlamydomonas reinhardtii Selected references

116.4 | 130 kDa

Chung et al. (2019) Distinct roles of Argonaute in the green alga Chlamydomonas reveal evolutionary conserved mode of miRNA-mediated gene expression. Sci Rep. 2019 Jul 31;9(1):11091. doi: 10.1038/s41598-019-47415-x.

Application example



Total proteins of Chlamydomonas reinhardtii saturated in 8M urea were separated on 15% SDS-PAGE and blotted for 1hour to 0.2 µm nitrocellulose at 100V using wet transfer system. Blots were blocked with 0.5% cold fish gelatin for 1hr at room temp with agitation. Blot was incubated in the primary antibody (anti-AGO3) at a dilution of 1:2500 for an hour at RT with agitation. The blots were washed with 3X 15min TBS-TT at RT with agitation. Blots as incubated in the secondary antibody (anti-rabbit, DyLight® 800 conjugated, AS12 2460 Agrisera) 1:5000 dilution for 30min at RT with agitation and washed 1X with TBSTT for 15min, 1X with TBST for 15min before scanning with the ODyssey IRD scanner

Courtesy Dr. Betty Chung, Department of Plant Sciences, University of Cambridge, UK