

This product is for research use only (not for diagnostic or therapeutic use)

contact: support@agrisera.com

Agrisera AB | Box 57 | SE-91121 Vännäs | Sweden | +46 (0)935 33 000 | www.agrisera.com

## Product no AS08 277

## Anti-APC | allophycocyanin alpha and beta

## **Product information**

Immunogen I native allophycocyanin alpha and beta purified from *Porphyridium cruentum* phycobilisomes

Host Rabbit

Clonality Polyclonal

Purity Serum

Format Lyophilized

Quantity 50 μl

**Reconstitution** For reconstitution add 50 μl of sterile water

Store lyophilized/reconstituted at -20°C; once reconstituted make aliquots to avoid repeated freeze-thaw cycles. Please remember to spin the tubes briefly prior to opening them to avoid any losses that might occur from material adhering to

the cap or sides of the tube

Additional information This product can be sold containing proclin if requested

## Application information

**Recommended dilution** 1:200-1:500 (IL), 1:1500 - 3000 (WB)

Expected | apparent 14-19 kDa

Confirmed reactivity Porphyridium cruentum, Synechocystis PCC 6803

Predicted reactivity Red Algae, Cyanobacteria

Species of your interest not listed? Contact us

Not reactive in No confirmed exceptions from predicted reactivity are currently known

Selected references Ge at al. (2017). Translating Divergent Environmental Stresses into a Common Proteome Response through the Histidine Kinase 33 (Hik33) in a Model Cyanobacterium. Mol Cell Proteomics. 2017 Jul;16(7):1258-1274. doi:

10.1074/mcp.M116.068080.

Gandini et al. (2017). The transporter SynPAM71 is located in the plasma membrane and thylakoids, and mediates manganese tolerance in Synechocystis PCC6803. New Phytol. 2017 Mar 20. doi: 10.1111/nph.14526.

Gunnelius et al. (2014). The omega subunit of the RNA polymerase core directs transcription efficiency in cyanobacteria. Nucleic Acids Res. 2014 Jan 29.

Hernandez-Prieto et al. (2011). The small CAB-like proteins of the cyanobacterium Synechocystis sp. PCC 6803: Their involvement in chlorophyll biogenesis for Photosystem II. Bioch.Bioph. Acta.

Gantt & Lipschultz (1974). Phycobilisome structure by immuno-electron microscopy. J. Phycology, Vol. 13:3, pages: 185-192. (immunolocalization)

Gantt & Lipschultz (1974). Phycobilisomes of Porphyridium cruentum: Pigment Analysis. Biochem. 13:2960. Gantt E & C Lipschultz (1977). Probing phycobilisome structure by immuno-electron microscopy. J Phycol. 13:18