

This product is for research use only (not for diagnostic or therapeutic use)

contact: support@agrisera.com

Agrisera AB | Box 57 | SE-91121 Vännäs | Sweden | +46 (0)935 33 000 | www.agrisera.com

Product no AS19 4322

## CAH6 | Carbonic anhydrase 6

## **Product information**

Immunogen Recombinant, full length CAH6 of Chlamydomonas reinhardtii, excised from the gel, UniProt: Q6S7R9

**Host** Rabbit

Clonality Polyclonal

Purity Serum

Format Lyophilized

Quantity 50 μl

**Reconstitution** For reconstitution add 50 μl, of sterile water

Storage Store lyophilized/reconstituted at -20°C; once reconstituted make aliquots to avoid repeated freeze-thaw cycles. Please remember to spin the tubes briefly prior to opening them to avoid any losses that might occur from material adhering to

the cap or sides of the tube.

## **Application information**

Recommended dilution 1:1000 (WB)

Expected | apparent

11 28 kDa

Confirmed reactivity Chlamydomonas reinhardtii

Predicted reactivity Gonium pectorale, Volvox carteri

Species of your interest not listed? Contact us

Not reactive in Lycopersicum esculentum, Neopyropia yezoensis (macroalga), Zea mays

Additional information Lately localization to flagella has been confirmed for CAH6 Mackinder et al. 2017.

Extraction method – harvest the cells and add SDS sample buffer.

Selected references Mackinder et al. 2017. A Spatial Interactome Reveals the Protein Organization of the Algal CO2-Concentrating

Mechanism. Cell. 2017 Sep 21;171(1):133-147.e14. doi: 10.1016/j.cell.2017.08.044.

Mitra et al. (2004). Identification of a new chloroplast carbonic anhydrase in Chlamydomonas reinhardtii. Plant Physiol.

2004 May;135(1):173-82.