

This product is for research use only (not for diagnostic or therapeutic use)

contact: support@agrisera.com

Agrisera AB | Box 57 | SE-91121 Vännäs | Sweden | +46 (0)935 33 000 | www.agrisera.com

### Product no AS09 530

# Anti-CBP20 | nuclear cap-binding protein subunit 2

### **Product information**

Immunogen KLH-conjugated peptide, derived with Arabidopsis thaliana CBP20 protein Q9xFD1, At5q44200

**Host** Rabbit

Clonality Polyclonal

**Purity** Immunogen affinity purified serum in PBS pH 7.4.

Format Lyophilized

Quantity 200 μg

**Reconstitution** For reconstitution add 200 μl of sterile water

Store lyophilized/reconstituted at -20°C; once reconstituted make aliquots to avoid repeated freeze-thaw cycles. Please remember to spin the tubes briefly prior to opening them to avoid any losses that might occur from material adhering to the cap or sides of the tube.

## Application information

Recommended dilution 1:1000 (WB)

Expected | apparent

29.6 | 30 kDa

Predicted reactivity

Glycne max, Hordeum vulgare, Lotus corniculatus, Nicotiana tabacum, Oryza sativa, Ricinus communis, Solanum lycopersicum, Solanum tuberosum, Zea mays

Species of your interest not listed? Contact us

Not reactive in No confirmed exceptions from predicted reactivity are currently known

Selected references

Raczynska et al. (2013). The SERRATE protein is involved in alternative splicing in Arabidopsis thaliana. Nucleic Acids Res. Oct 16.

#### **Application example**



25 µg of total protein extratcs from 10 days old seedlings of Arabidopsis thaliana (wild type and a CBP20 mutant) were separated on 12.5 % SDS-PAGE and blotted 1h to PVDF (tank blotting). Blots were blocked with Roti-block over night at 4°C agitation. Blot was incubated in the primary antibody at a dilution of 1:1 000 for 1h at RT with agitation. The antibody solution was decanted and the blot was rinsed briefly twice, then washed once for 15 min and 3 times for 5 min in TBS-T at RT with agitation. Blot was incubated in secondary antibody (anti-rabbit IgG horse radish peroxidase conjugated, from Agrisera AS09 602) diluted to 1:10 000 in for 1h at RT with agitation. The blot was washed as above and developed for 5 min with ECL according to the manufacturers instructions. Exposure time was 30 min.

Courtesy of Dr. Sascha Laubinger, ZMBP, Germany