

This product is for research use only (not for diagnostic or therapeutic use)

contact: support@agrisera.com

Agrisera AB | Box 57 | SE-91121 Vännäs | Sweden | +46 (0)935 33 000 | www.agrisera.com

Product no AS16 3114

Anti-Fdx2 | Ferredoxin 2

Product information

Immunogen KLH-conjugated synthetic peptide chosen from a C-terminal of Chlamydomonas reinhardtii ferredoxin 2 Q2HZ25

Host Rabbit

Clonality Polyclonal

Purity Serum

Format Lyophilized

Quantity 50 ul

Reconstitution For reconstitution add 50 μl of sterile water

Storage Store lyophilized/reconstituted at -20°C; once reconstituted make aliquots to avoid repeated freeze-thaw cycles. Please remember to spin the tubes briefly prior to opening them to avoid any losses that might occur from material adhering to

the cap or sides of the tube.

Application information

Recommended dilution 1:1000 (WB)

Expected | apparent 13.4 kDa

MW

Confirmed reactivity Chlamydomonas reinhardtii

Predicted reactivity Chlamydomonas reinhardtii

Not reactive in No confirmed exceptions from predicted reactivity are currently known

Additional information Antibody is recognizing 10 ng of recombinant ferredoxin 2.

Load over 30 ug/well is recommended for soluble Chlamydomonas reinhardtii cell extract.

Terrauchi et al. (2009). Pattern of expression and substrate specificity of chloroplast ferredoxins from Chlamydomonas Selected references

reinhardtii. J Biol Chem 284 (38):25867-25878.