

This product is for research use only (not for diagnostic or therapeutic use)

contact: support@agrisera.com

Agrisera AB | Box 57 | SE-91121 Vännäs | Sweden | +46 (0)935 33 000 | www.agrisera.com


Product no AS13 2651 Anti-FtsZ2 | Plant cell division protein ftsZ2

Product information

| Immunogen | Recombinant part of Arabidopsis thalina FtsZ2 conserved in FtsZ2-1 UniProt: <u>082533.</u> TAIR: <u>At2g36250</u> , FtsZ2-2 UniProt: <u>Q9LXJ0</u> , TAIR: <u>AT3G52750</u> | |
|-------------------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|--|
| Host | Rabbit | |
| Clonality | Polyclonal | |
| Purity | Serum | |
| Format | Lyophilized | |
| Quantity | 50 μl | |
| Reconstitution | For reconstitution add 50 μ l of sterile water | |
| Storage | Store lyophilized/reconstituted at -20°C; once reconstituted make aliquots to avoid repeated freeze-thaw cycles. Please remember to spin the tubes briefly prior to opening them to avoid any losses that might occur from material adhering to the cap or sides of the tube. | |
| Additional information | Antibody available in 3 various pack sizes: 50, 100 and 150 μ l - please, inquire | |
| Application information | | |
| Recommended dilution | 1 · 2000 (WB) | |

| Recommended dilution | 1 : 2000 (WB) |
|---------------------------|-----------------------------------------------------------------------|
| Expected apparent MW | 50.7 kDa |
| Confirmed reactivity | Arabidopsis thaliana, Hordeum vulgare |
| Predicted reactivity | Glycne max, Jatropha manihot, Ricinus communis |
| | Species of your interest not listed? Contact us |
| Not reactive in | No confirmed exceptions from predicted reactivity are currently known |

Application example



2 μg of recombinant FtsZ2 fragment(~10 kDa) (1) and total protein extracted with PEB <u>AS08 300</u> from leaf tissue *Arabidopsis thaliana* (2) were separated on 4-12% NuPage (Invitrogen) LDS-PAGE and blotted 1h to PVDF. Filters were blocked 1h with 2% low-fat milk powder in TBS-T (0.1% TWEEN 20) and probed with FtsZ2 antibody (AS13 2651, 1:2000, 1h) and secondary anti-rabbit (1:25 000, 1 h) antibody (HRP conjugated, secondary antibody <u>AS09 602</u> in TBS-T containing 2% low fat milk powder. Antibody incubations were followed by washings in TBS-T (15, +5, +5, +5 min). All steps were performed at RT with agitation. The signal was developed for 5 min with ECL according to the manufacturer's instructions. Exposure time was 2 min.