

This product is for research use only (not for diagnostic or therapeutic use)

contact: support@agrisera.com

Agrisera AB | Box 57 | SE-91121 Vännäs | Sweden | +46 (0)935 33 000 | www.agrisera.com

Product no AS20 4443

Anti-GFP | Green Fluorescence Protein

Product information

Immunogen Recombinant GFP protein derived from Aequorea victoria, UniProt: P42212

Host Rabbit

Clonality Polyclonal

Purity Immunogen affinity purified serum in PBS pH 7.4.

Format Lyophilized

Quantity 50 μg

Reconstitution For reconstitution add 50 μl of sterile water

Storage Store lyophilized/reconstituted at -20°C; once reconstituted make aliquots to avoid repeated freeze-thaw cycles. Please remember to spin the tubes briefly prior to opening them to avoid any losses that might occur from material adhering to

the cap or sides of the tube.

Application information

Recommended dilution 1: 5000 - 1 : 10 000 (WB)

Confirmed reactivity Recombinant GFP overexpressed in *E.coli* from *Arabidopsis thaliana, Chlamydomonas reinhardii, Nicotiana tabacum*

Predicted reactivity mClover3

Selected references Dieren et al. (2024). Analysis of abiotic and biotic stress-induced Ca2+ transients in the crop species Solanum

tuberosum. Sci Rep . 2024 Nov 11;14(1):27625. doi: 10.1038/s41598-024-79134-3.



19.4 μg/well of total protein extracted freshly from tobacco with buffer (50mM Tris pH 7.5, 150 mM NaCl, 1 mM EDTA, 10% Glicerol, 1 mM DTT, 1X Pefablock) and denatured with SDS 2% and 0.1% Bromophenol blue (Laemmli buffer) at 95°C for 5 min. Separated on 12 % SDS-PAGE and blotted 2h to PVDF/nitrocellulose, using wet transfer. Blot was blocked with 5% milk for 1h/RT with agitation. Blot was incubated in the primary antibody at a dilution of 1: 10 000 in TBS-T ON/4°C with agitation. The antibody solution was decanted and the blot was rinsed briefly twice, then washed 3 times for 10 min in TBS-T at RT with agitation. Blot was incubated in Agrisera matching secondary antibody (anti-rabbit IgG horse radish peroxidase conjugated, <u>AS09 602</u>) diluted to 1:20 000 in for 1h/RT with agitation. The blot was washed as above and developed for 2 min with <u>Agrisera ECLSuperBright</u>. Exposure time was 60 seconds.


Courtesy of Delfina Gagliardi Instituto de Agrobiotecnología del Litoral (IAL), Argentina



This product is **for research use only** (not for diagnostic or therapeutic use)

contact: support@agrisera.com

Agrisera AB | Box 57 | SE-91121 Vännäs | Sweden | +46 (0)935 33 000 | www.agrisera.com



Super-resolution image of immunofluorescent localization of GFP-fused protein in *Nicotiana tabacum* (tobacco) leaf tissue using Agrisera anti-GFP antibodies (AS20 4443) and donkey anti-rabbit IgG DyLight 650 conjugated secondary antibodies (Green). Calcofluor White (Blue) and Tissue Autofluorescence (Yellow),

Method:

Fixation: 2% paraformaldehyde + 0.5% glutaraldehyde, dehydration in ethanol, embedment and UV polymerization in HM20 resin

Cell wall digestion: no

Membrane permeabilization: DMSO-IGEPA

Air drying; No Antigen retrieval: No

Blocking buffer: 4% w/v non-fat milk powder in 1X Phosphate Buffered Saline buffer, prepared freshly, filtered with 0.2 µm filer

Washing buffer: 1X Phosphate Buffered Saline buffer Primary antibody dilution and incubation time: 1:100; 1 hour at room temperature

Secondary antibody dilution and incubation time and supplier: 1:200 Donkey anti-rabbit IgG (H&L) DyLight®650 (Agrisera, AS12 2254); 1 hour at recommendation

room temperature Co-staining of the nucleus (DAPI): No

Cell wall and nucleus staining: Calcofluor White

All material was freshly prepared. All incubation volumes were kept the same: 200 µL per 22x22 mm cover glass/sample

Courtesy of Dr. Kirk Czymmek, Donald Danforth Plant Science Centre, USA.