

This product is for research use only (not for diagnostic or therapeutic use)

contact: support@agrisera.com

Agrisera AB | Box 57 | SE-91121 Vännäs | Sweden | +46 (0)935 33 000 | www.agrisera.com

## Product no AS15 2914

## Anti-GST conjugated to SureLight® R-Phycoerythrin

## **Product information**

**Immunogen** Glutathione-S-transferase from *S. japonicum* 

Host Goat

Clonality Polyclonal

Purity Immunogen affinity purified using solid phase GST in 100 mM Sodium Phosphate (pH 7.4), 100 mM NaCl, 50 mM

Sucrose, 1.5% BSA and 0.01% Thimerosal.

Format Lyophilized

Quantity 100 μg

**Reconstitution** For reconstitution add 100 μl of sterile water

Storage Storag

add an equal volume of glycerol and Store at -20°C.

Additional information Antibody is provided in: 100 mM Sodium Phosphate (pH 7,4), 100 mM Sodium Chloride, 50 mM Sucrose, 0,1% BSA

and 2 mM Sodium Azide

## **Application information**

**Recommended dilution** The optimal working dilution should be determined by the investigator