

This product is for research use only (not for diagnostic or therapeutic use)

contact: support@agrisera.com

Agrisera AB | Box 57 | SE-91121 Vännäs | Sweden | +46 (0)935 33 000 | www.agrisera.com

Product no AS99 001

Anti-KLH | Keyhole limpet hemocyanin

Product information

Immunogen Purified keyhole limpet hemocyanin (<u>KLH</u>), whole molecule,

Host Rabbit

Clonality Polyclonal

Purity Immunogen affinity purified serum in PBS pH 7.4.

Format Lyophilized

Quantity 50 μg

Reconstitution For reconstitution add 50 μl of sterile water

Storage Store lyophilized/reconstituted at -20 °C; once reconstituted make aliquots to avoid repeated freeze-thaw cycles. Please

remember to spin the tubes briefly prior to opening them to avoid any losses that might occur from material adhering to

the cap or sides of the tube.

Additional information Protein present in plant vascular tissue (xylem and vascular cambium) is detected by anti-KLH antibodies (Höglund et al. 2002) which might lead to false results in IL when using anti-peptide antibodies generated to KLH-conjugated

ai. 2002) which might lead to false results in it when using anti-peptide antibodies generated to KLH-peptide.

Application information

Recommended dilution 1: 10 000 (ELISA), 1: 10 000 (WB), 1: 1000 (IL)

Expected | apparent ca. 400 kDa/subunit

Confirmed reactivity Megathura crenulata - most commonly used carrier protein

Not reactive in No confirmed exceptions from predicted reactivity are currently known

Additional information Antibody can be used as a negative control to determine if observed signal is generated by anti-KLH or anti-peptide

antibodies, Due to its large size KLH protein will be very difficult to separate on SDS-PAGE

Selected references Geadkaew et al. (2014). Bi-functionality of Opisthorchis viverrini aquaporins. doi:10.1016/j.biochi.2014.11.013.

Höglund et al. (2002). An Antigen Expressed During Plant Vascular Development Crossreacts with Antibodies Towards

KLH (Keyhole Limpet Hemocyanin). J of Histochem & Cytochem. 50:999-1003.