

This product is for research use only (not for diagnostic or therapeutic use)

contact: support@agrisera.com

Agrisera AB | Box 57 | SE-91121 Vännäs | Sweden | +46 (0)935 33 000 | www.agrisera.com

Product no AS18 4197-1ml

Anti-Pectic polysaccharide, beta-1,4-galactan (monoclonal, clone LM5)

Product information

Immunogen Pectic polysaccharide (beta 1,4-galactan)

Host Rat

Clonality Monoclonal

Subclass/isotype IgG

Purity Cell culture supernatant.

Format Liquid

Quantity 1 ml

Storage

Store at +4°C (short term) and at -20°C (long term). Make aliquots to avoid repeated freeze-thaw cycles. Please remember to spin the tubes briefly prior to opening them to avoid any losses that might occur from any material

adhering to the cap or sides of the tube.

Additional information

Contains 0.05% Sodium Azide.

No cross-reactivity with (1-3)-beta-D-galactans or (1-6)-beta-D-galactans.

It recognizes a linear tetrasaccharide in (1-4)-beta-D-galactans.

In ELISA (competitive inhibition), antibody is binding to: (1-4)-beta-D-galactan was inhibited (50%) by 58 µg/ml (1-4)-beta-D-galactotetraose and by 0.7 $\mu g/ml$ lupin (1-4)-beta-D-galactan.

Application information


Recommended dilution 1:10 (ELISA, IF)

Confirmed reactivity | Higher plants, ferns and mosses

Not reactive in No confirmed exceptions from predicted reactivity are currently known

Selected references

Andersen et al. (2016). Characterization of the LM5 pectic galactan epitope with synthetic analogues of ß-1,4-d-galactotetraose. Carbohydr Res. 2016 Dec ;436:36-40.doi: 10.1016/j.carres.2016.10.012. Jones et al. (1997). Development and validation of an in vitro model system to study peripheral sensory neuron development and injury. Sci Rep. 2018 Oct 29;8(1):15961. doi: 10.1038/s41598-018-34280-3.



Samples:

- A Lignin autofluorescence
- B LM5 staining of wood
- C Composite
- D LM5 antibody staning + brightfield Objective: 20x,100 µm scalebar

Type of material: Pine (Pinus sylvestris) reaction wood section in paraffin cut with microtome

Fixation: none Hydrophilization: no Cell wall digestion: no

Membrane permeabilization: none Antigen retrieval: no

Blocking buffer: 5% milk in PBS Washing buffer: 5% milk in PBS-T



This product is for research use only (not for diagnostic or therapeutic use)

contact: support@agrisera.com

Agrisera AB | Box 57 | SE-91121 Vännäs | Sweden | +46 (0)935 33 000 | www.agrisera.com

Primary antibody dilution and incubation time: 1:20 incubation 90 minutes, room temperature

Secondary antibody dilution and incubation time and supplier: 1:100 Goat anti-Rat IgG (H+L) Cross-Adsorbed Secondary Antibody, Alexa Fluor™ 594 (invitrogen: A-11007), incubation for 60 minutes in room temperature.

Co-staining of the nucleus (DAPI): no Cell wall and nucleus staining: none

Courtesy of Dr. Jan Łyczakowski, Jagiellonian University, Kraków, Poland