

This product is for research use only (not for diagnostic or therapeutic use)

contact: support@agrisera.com

Agrisera AB | Box 57 | SE-91121 Vännäs | Sweden | +46 (0)935 33 000 | www.agrisera.com

Product no AS22 4730

Anti-PsaD | PSI-D subunit of photosystem I (algal)

Product information

Immunogen PsaD of Chlamydomonas reinhardtii, UniProt: Q39615

Host Rabbit

Clonality Polyclonal

Purity Serum

Format Lyophilized

Quantity 50 μl

Reconstitution For reconstitution add 50 μl of sterile water

Storage Store lyophilized/reconstituted at -20 °C; once reconstituted make aliquots to avoid repeated freeze-thaw cycles,

Please, remember to spin tubes briefly prior to opening them to avoid any losses that might occur from lyophilized

material adhering to the cap or sides of the tubes

Application information

Recommended dilution 1: 5000 - 1 : 10 000 (WB)

Expected | apparent 21.3 kDa

MW

Confirmed reactivity | Chlamydomonas reinhardtii

Additional information Can be provided with ProClin, if requested