## 

This product is for research use only (not for diagnostic or therapeutic use)

Agrisera AB | Box 57 | SE-91121 Vännäs | Sweden | +46 (0)935 33 000 | www.agrisera.com

## Product no AS22 4732 Anti-PsaF | PSI-F subunit of photosystem I (algal)

## **Product information**

 Immunogen
 PsaF of Chlamydomonas reinhardtii, UniProt: <u>A8J4S1</u>

 Host
 Rabbit

 Clonality
 Polyclonal

 Purity
 Serum

 Format
 Lopphilized

 Quantity
 So µl

 Reconstitution
 For reconstitution add 50 µl of sterile water

 Storage
 Store Iyophilized/reconstituted at -20°C; once reconstituted make aliquots to avoid repeated freeze-thaw cycles, Please, remember to spin tubes briefly prior to opening them to avoid any losses that might occur from lyophilized material adhering to the cap or sides of the tubes

## **Application information**

Recommended dilution	1 : 5 000 (WB)
Expected   apparent MW	24 kDa
Confirmed reactivity	Chlamydomonas reinhardtii
Additional information	Can be provided with ProClin, if requested

1/1