

This product is for research use only (not for diagnostic or therapeutic use)

contact: support@agrisera.com

Agrisera AB | Box 57 | SE-91121 Vännäs | Sweden | +46 (0)935 33 000 | www.agrisera.com

Product no AS21 4554

Anti-Rad22 | DNA repair and recombination protein rad22 (Saccharomyces pombe)

Product information

Immunogen Purified, full length, recombinant Rad22 protein from Saccharomyces cerevisiae, UniProt: P36592 overexpressed in E.

Host Rabbit

Clonality Polyclonal

Purity Immunogen affinity purified serum, in PBS. Contains 50 % glycerol, filter sterilized.

Format Liquid

Quantity 50 ug

Storage Store at -20 °C; make aliquots to avoid repeated freeze-thaw cycles. Please remember to spin the tubes briefly prior to opening them to avoid any losses that might occur from material adhering to the cap or sides of the tube.

Application information

Recommended dilution 1:2000 - 1:5000 (WB)


Expected | apparent 52 | 58 kDa

Confirmed reactivity Saccharomyces pombe


Predicted reactivity Species of your interest not listed? Contact us

Selected references Lehmann (1996). Molecular biology of DNA repair in the fission yeast Schizosaccharomyces pombe. Mutat Res. 1996

Aug 8;363(3):147-61. doi: 10.1016/0921-8777(96)00017-1. PMID: 8765156.



Whole cell extract of Saccharomyces pombe was separated on SDS-PAGE and blotted to a membrane. Primary antibody was incubated at 1: 2000, followed by washes and incubation with a secondary goat anti-rabbit IgG HRP conjugated antibodies, used at 1: 10 000 1h/RT. Reaction was developed using chemiluminescence following manufacture's recommendations.



Whole cell extract of Saccharomyces pombe wild-type (1) and Saccharomyces pombe Rad22 deletion mutant (2) was separated on a 12.5 % SDS-PAGE and blotted to a membrane using wet transfer. Primary antibody was incubated at 1: 2000, followed by washes and incubation with a secondary goat anti-rabbit IgG HRP conjugated antibodies, used at 1: 10 000 1h/RT. Reaction was developed using chemiluminescence following manufacture's recommendations.