

This product is for research use only (not for diagnostic or therapeutic use)

contact: support@agrisera.com

Agrisera AB | Box 57 | SE-91121 Vännäs | Sweden | +46 (0)935 33 000 | www.agrisera.com

Product no AS21 4687 Anti-TAP epitope tag

Product information

Immunogen	TAP epitope tag, sequence: CSSGALDYDIPTTASENLYFQ, derived from the C-terminus of the TAP-tag construct after TEV cleavage,
Host	Mouse
Clonality	Monoclonal
Purity	Affinity purified in PBS pH 7.4. Contains 0.02 % sodium azide. Contains 50 % glycerol.
Format	Liquid
Quantity	50 µg
Storage	Store at -20°C, Make aliquots to avoid repeated freeze-thaw cycles. Please remember to spin the tubes briefly prior to opening them to avoid any losses that might occur from material adhering to the cap or sides of the tube.
ation information	

Application information Recommended dilution 1: 500 - 1: 5000 (WB) Expected | apparent MW Depends upon protein fusion partner Confirmed reactivity TAP epitope tag Selected references To be added when available, antibodies released in October 2021.