

Product no **AS21 4687**

TAP epitope tag

Product information

| | |
|------------------|---|
| Immunogen | TAP epitope tag, sequence: CSSGALDYDIPTTASENLYFQ, derived from the C-terminus of the TAP-tag construct after TEV cleavage, |
| Host | Mouse |
| Clonality | Monoclonal |
| Purity | Affinity purified in PBS pH 7.4. Contains 0.02 % sodium azide. Contains 50 % glycerol. |
| Format | Liquid |
| Quantity | 50 µg |
| Storage | Store at -20 °C, Make aliquots to avoid repeated freeze-thaw cycles. Please remember to spin the tubes briefly prior to opening them to avoid any losses that might occur from material adhering to the cap or sides of the tube. |

Application information

| | |
|-------------------------------|--|
| Recommended dilution | 1: 500 - 1: 5000 (WB) |
| Expected apparent MW | Depends upon protein fusion partner |
| Confirmed reactivity | TAP epitope tag |
| Selected references | To be added when available, antibodies released in October 2021. |